

Accompanying the Recovery Journey

Create a Safer Environment

YUWELL 7A-23B ELECTRIC SUCTION MACHINE

SKU: PRO195

www.procaredmedfurn.com

- **Required Input Power:** 180 VA at 50 Hz, AC 220 V
- **Pumping Capability:** 40 LPM of pumping capacity
- **Negative Power:** 90 kPa is the maximum negative pressure.
- **Noise Level Of the machine:** Noise Level at Operation: < 65 dB
- **Jar Capacity:** 2.5 L x 2 Principal Attributes:

KEY FEATURE


- THESE ELECTRIC SUCTION MACHINES OFFER MAINTENANCE AND OIL-FREE SUCTION PUMPS FOR THE BEST SUCTION.
- COMES WITH USER-FRIENDLY AND SMOOTH CONTROLS
- THE OUTER BODY OF THIS ELECTRIC SUCTION MACHINE IS MADE FROM A PREMIUM META WITH POWDER COATING
- FOR THE BEST EASIEST AND MOST ACCURATE OPERATION, USE FOOTSWITCH, TUBING SYSTEM, AND EXTRA FILTERS.
- THESE MACHINE ARE PERFECT FOR HOSPITALS, AND EMERGENCIES WHERE EFFECTIVE SUCTION IS NEEDED.



CE

Procure

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaredmedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procure