

Accompanying the Recovery Journey

Create a Safer Environment

YUWELL 7E-B5

SUCTION MACHINE

SKU: PRO200

www.procaredmedfurn.com

Power Input: 90 VA, AC 220 V at 50 Hz

Pumping Rate: 22 L/min

Negative pressure – 75 kPa.

Noise Level at Operation: < 45 dB

Capacity of Jar: 1000 mL

KEY FEATURE

VERSATILE USE:

7E-B5 IS A PERFECT MACHINE YOU CAN USE THIS MACHINE IN DIFFERENT AREAS. YOU CAN USE IN-HOME CARE, CLINICS, OPERATING ROOMS, AND IN EMERGENCIES.

INNOVATIVE DESIGN:

ITS DESIGN MAKE IT MORE COMPACT AND SIMPLE FOR TRANSPORTATION, THE MACHINE COMES WITH A STRONG AND LIGHTWEIGHT PLASTIC BODY

HIGH-PERFORMANCE PUMP:

7E-B5 OFFERS A MASSIVE 22 LPM PUMPING CAPACITY, WITH AN OIL-FREE, HIGH-PRESSURE SUCTION PUMP WHICH IS MAINTENANCE-FREE.

SILENT OPERATION:

TO MAINTAIN SILENCE IN THE HOSPITAL AREA WHICH IS ESSENTIAL FOR PATIENT COMFORT, WITH AN OPERATING NOISE LEVEL OF LESS THAN 45 DB

ADJUSTABLE PRESSURE:

THE NEGATIVE PRESSURE OF THIS SUCTION MACHINE CAN BE CONTROLLED USING A REGULATOR UP TO 75 KPA, THE BEST CONTROL ON YOUR PROCEDURE FOR ACCURATE SUCTION.

ENHANCED SAFETY:

TO MAINTAIN STRICT HYGIENE DURING WORK, A TUBING KIT, EXTRA BACTERIAL FILTERS, AND A 1000 ML POLYCARBONATE JAR WITH AN ADVANCED OVERFLOW PREVENTION SYSTEM ARE PROVIDED

SIMPLE SETUP:

THE DETACHABLE POWER CORD IS PROVIDED IT MAKES SETUP AND TRANSPORTATION EASIER AND USER-FRIENDLY.



Jar




Grab Handle



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaredmedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare