



# Accompanying the Recovery Journey

Create a Safer Environment

## Z TROLLEY

SKU: PRO104


[www.procaremedfurn.com](http://www.procaremedfurn.com)

- Length: 500 mm
- Width: 400 mm
- Height: 750 mm

## KEY FEATURE

- **UNIQUE DESIGN:**
- THIS MONITOR TROLLEY COMES IN A UNIQUE Z SHAPE DESIGN WHICH IMPROVES THE APPEARANCE OF THE HOSPITAL WARDS.
- **STRONG CONSTRUCTION:**
- THE Z-SHAPE MONITOR TROLLEY IS MADE FROM PREMIUM CRCA STEEL THAT PROVIDES DURABILITY, AND LONG LIFE TO THE PRODUCT. THIS MONITOR TROLLEY IS PRETREATED AND EPOXY POWDER-COATED MAKING THIS TROLLEY CORROSION-RESISTANT AND EASY TO CLEAN TO MAINTAIN THE HYGIENE OF THE WARDS.
- **MULTIPLE SHELVES:**
- THIS MONITOR TROLLEY COMES WITH THREE SHELVES AND OFFERS A SPACE TO KEEP MONITORS, MEDICINES, OR MEDICAL SUPPLIES.
- **WHEELS:**
- THE TROLLEY IS EQUIPPED WITH 50 MM PREMIUM WHEELS FOR EASY MOVEMENT, AND TO CARRY THE TROLLEY EASILY BETWEEN THE WARDS.
- 

Powder  
Coated



Wheel



CE

**Procare**

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare