



# Accompanying the Recovery Journey

Create a Safer Environment

## BEDSIDE DRUG TROLLEY WITH SS TOP

SKU: PRO119B

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- Length: 580 mm
- Width: 360 mm
- Height: 660 mm


### KEY FEATURE

- **STAINLESS STEEL TOP:**
- THIS DRUG TROLLEY COMES WITH A STAINLESS STEEL TOP SURFACE, SS SURFACES ARE EASY TO CLEAN AND MAINTAIN HYGIENIC CONDITIONS IN THE WARDS.
- **STRONG DESIGN:**
- THE DRUG TROLLEY IS PRETREATED AND EPOXY POWDER COATED. WHICH MAKES IT CORROSION-RESISTANT, DURABLE, AND USEFUL IN DAILY HARSH CONDITIONS. ALSO AVAILABLE IN FULLY STAINLESS STEEL CONSTRUCTION.
- **MULTI SHELVES:**
- YOU CAN CUSTOMIZE YOUR TROLLEY AS YOU WANT, THE TROLLEY COMES WITH THREE SHELVES, AND YOU CAN CUSTOMIZE IT TO FOUR SHELVES
- **MOBILITY:**
- THIS BEDSIDE DRUG TROLLEY OFFERS SMOOTH AND STABLE MOVEMENT WITH ITS 50 MM DIAGONAL WHEELS.
- 

Powder  
Coated



Stainless  
Steel



Grab  
Handle



Wheel



CE

## Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare