



Accompanying the Recovery Journey

Create a Safer Environment

DOUBLE OXYGEN CYLINDER TROLLEY www.procaremedfurn.com

- Length – 660 mm
- Width – 610 mm
- Height – 1400 mm

KEY FEATURE

PREMIUM CONSTRUCTION:
THIS OXYGEN CYLINDER TROLLEY IS MADE FROM PREMIUM CRCA STEEL, AND COATED WITH EPOXY POWDER. ITS PREMIUM CONSTRUCTION MAKES THIS TROLLEY MORE DURABLE, AND RUST RESISTANT AND OFFERS LONG LIFE TO THE PRODUCT

SAFETY CHAINS:
THE SAFETY CHAINS ARE PROVIDED AT THE BOTTOM AND TOP OF THE TROLLEY TO HOLD CYLINDERS, THIS SAFETY CHAINS PREVENT THE FALL OF THE CYLINDER FROM THE TROLLEY DURING TRANSPORTATION.

PREMIUM WHEELS:
PREMIUM WHEELS ARE PROVIDED FOR SMOOTH AND EASY TRANSPORTATION ON VARIOUS FLOORS. THIS WHEEL HELPS TO PUSH THE TROLLEY WITHOUT EXTRA EFFORT.

EASY TO LIFT :
ITS FRONT SIDE EDGE HELPS TO LIFT THE HEAVY CYLINDER EASILY.

Grab Handle



Safety Chain



Wheel



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare