



# Accompanying the Recovery Journey

Create a Safer Environment

HORIZONTAL  
AUTOCLAVE – FULLY  
AUTOMATIC  
[www.procaremedfurn.com](http://www.procaremedfurn.com)

## KEY FEATURE

- CAPACITY OF THE AUTOCLAVE – THIS HIGH-PRESSURE AUTOCLAVE COMES WITH A 1000 LITER CAPACITY, IT'S IDEAL FOR HIGH-VOLUME STERILIZATION
- CONSTRUCTION: THE CONTACT PARTS ARE MADE FROM STAINLESS STEEL 304 Q FOR THE BEST DURABILITY AND TO MAKE IT CORROSION RESISTANT. AND NON-CONTACT PARTS MADE FROM MS.
- DOOR DESIGN: THIS HORIZONTAL AUTOCLAVE COMES WITH A SINGLE-DOOR DESIGN FOR EASY AND SAFE OPERATION.
- CONTROL PANEL: THE AUTOCLAVE IS EQUIPPED WITH A TEMPERATURE AND TIME CONTROLLER FOR ACCURATE OPERATION AND TEMPERATURE CONTROL
- OPERATING TEMPERATURE: THE TEMPERATURE WITH 121°C WITH WORKING PRESSURE OF 1.2 KG / CM2
- HEATING SYSTEM: THE AUTOCLAVE IS EQUIPPED WITH WATER IMMERSION HEATER ( 6 X 6 )
- STEAM CAPACITY: THE STEAM GENERATION CAPACITY OF THE AUTOCLAVE IS 36 KW WITH A TOTAL CAPACITY OF 42 KW.
- VACUUM PUMP AND CONDENSER: THE VACUUM PUMP AND CONDENSER HELP TO ENHANCE THE PERFORMANCE OF THE AUTOCLAVE FOR AN EFFICIENT STERILIZATION
- SWITCHGEAR: THE SWITCHGEAR OF THE AUTOCLAVE IS POWERED BY SCHNEIDER ELECTRICALS.
- TROLLEY FOR EFFICIENCY: THE LOADING AND UNLOADING TROLLEY HELPS TO LOAD AND UNLOAD EQUIPMENT EASILY FROM THE AUTOCLAVE MACHINE
- MANUAL OPERATION: YOU CAN CONTROL THIS AUTOCLAVE MACHINE MANUALLY USING INDIVIDUAL SWITCHES, IN CASE OF AUTOMATION FAILURE, TO PROVIDE



**Procare**

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

**Procare**