



PRO-ME-101

www.procaremedfurn.com

- Intensity: 25,000 Lux
- No. of LED's: 5
- Spot Diameter: 100-150 mm
- Power Supply: 220V/50Hz AC
- Length: 1.6mtr
- Power Consumption: 4 watt
- Colour Temp.: 3800 – 4800k
- LED Average life: >50000 Hrs.
- Focusing: Fix
- Battery Backup (Optional)

Accompanying the Recovery Journey

Create a Safer Environment



CE

ProCare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

ProCare