



Accompanying the Recovery Journey

RADIANT WARMER
WITH RESUSCITATION
UNIT
www.procaremedfurn.com

Create a Safer Environment

KEY FEATURE

- Introducing our state-of-the-art Radiant Warmer With Resuscitation Unit, a groundbreaking solution designed to provide unparalleled care for newborns in medical settings. This innovative system seamlessly combines the functionalities of a radiant warmer with essential features of a resuscitation unit, ensuring superior support during critical stages of infancy.
- Equipped with a fixed baby cradle and a premium-quality mattress, our warmer prioritizes comfort and stability for newborns. Driven by an advanced micro-processor-based servo Controller, it offers precise temperature control across various modes, including Servo skin, AIR, Pre-Warm, and Manual, delivering tailored thermal management to meet individual patient needs.
- Featuring dual digital displays for real-time monitoring of both baby and set temperatures, our warmer empowers healthcare providers with accurate environmental control. Intuitive soft-touch buttons and customizable temperature units (°F and °C) enhance usability and adaptability, ensuring seamless integration into clinical workflows.
- In addition to essential functionalities such as a Baby Feed Reminder, Monitor and IV stand, Swiveling top for X-ray and Easy Access, LED Illumination for observation, and safety alarms with high-temperature cut-off, our warmer boasts a versatile resuscitation unit for emergency interventions.
- Constructed with clear acrylic walls on all sides and collapsible front and side panels, the resuscitation unit offers exceptional visibility and accessibility during critical procedures. The baby tray, accompanied by a high-quality mattress, guarantees optimal comfort for newborns, while the clear acrylic base facilitates under-surface phototherapy.
- With three convenient drawers at the base for user convenience, our Radiant Warmer With Resuscitation Unit is a comprehensive solution for neonatal care in hospitals and medical facilities. Elevate the well-being and safety of newborns with our innovative and adaptable system.
-



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY