



BLOOD COLLECTION CHAIR WITH FOLDING FOOT REST

www.procaremedfurn.com

SIZE/DIMENSION :- 39LX20WX50H INCH

FRAME MATERIAL :- MS

USAGE/ APPLICATION :- HOSPITAL & CLINIC

COLOR :- CUSTOMIZED

Accompanying
the Recovery Journey

KEY FEATURE

THE DENTAL CHAIR CLASSIC IS A HIGH-PERFORMANCE AND ERGONOMICALLY DESIGNED DENTAL TREATMENT UNIT, OFFERING SEAMLESS FUNCTIONALITY AND ENHANCED PATIENT COMFORT. IT FEATURES ELECTRICAL OPERATION WITH TWO FIXED POSITIONS—ZERO AND SUPINE, ENSURING EFFORTLESS POSITIONING FOR DENTAL PROCEDURES. THE CHAIR'S SOFT, SEAMLESS UPHOLSTERY PROVIDES SUPERIOR PATIENT COMFORT, WHILE THE ABS-MOULDED ARMREST ADDS CONVENIENCE AND DURABILITY. EQUIPPED WITH A CERAMIC SPITTOON AND TUMBLER, IT INCLUDES AN AUTO-CLOSE WATER CUP FILLER AND SPITTOON FLUSH, PROMOTING HYGIENE AND EASY MAINTENANCE. THE 4-LED SENSOR LAMP WITH ADJUSTABLE ILLUMINATION FROM 7,000 TO 32,000 LUX ENSURES EXCELLENT VISIBILITY FOR PRECISION TREATMENTS.

FOR EFFICIENCY, THE CHAIR FEATURES A MULTIPURPOSE UTILITY TRAY DESIGNED TO HOLD ESSENTIAL TOOLS LIKE A SCALER, GLASS, AND NAPKIN, ALONG WITH A LARGE STAINLESS STEEL INSTRUMENT TRAY WITH A PNEUMATIC LOCKING SYSTEM FOR SECURE TOOL PLACEMENT. THE DELIVERY UNIT INCLUDES PROVISIONS FOR TWO AIROTORS, ONE MICROMOTOR, AND ONE THREE-WAY SYRINGE, WHILE THE ASSISTANT UNIT IS EQUIPPED WITH A THREE-WAY SYRINGE AND LOW SUCTION FOR SEAMLESS WORKFLOW. THE BUILT-IN X-RAY VIEWER WITH WHITE LIGHT AIDS IN CLEAR DIAGNOSIS. ADDITIONALLY, A MULTIFUNCTIONAL FOOT CONTROL AND DOUBLE-ARTICULATED HEADREST MAKE IT HIGHLY ADAPTABLE FOR PEDIATRIC AND HANDICAPPED PATIENTS.

WITH ITS MODERN FEATURES, ERGONOMIC DESIGN, AND USER-FRIENDLY CONTROLS, THE DENTAL CHAIR CLASSIC IS AN IDEAL CHOICE FOR DENTAL CLINICS LOOKING TO ENHANCE EFFICIENCY, PRECISION, AND PATIENT CARE.

Create a Safer Environment

Premium Quality
Cushion



Hand rest



Foot Rest



CE

Procare

SEC NO. 7 , PLOT NO. 146A, PCNTDA MIDC
BHOSARI, PUNE - 411039, MAHARASHTRA
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare