



DENTAL CHAIR AQUA

www.procaremedfurn.com

- ▶ ELECTRICALLY OPERABLE DENTAL TREATMENT UNIT WITH 2 FORED PROGRAMS: ZERO & SUPINE POSITION
- ▶ SEAMLESS & SOFT UPHOLSTERY FOR PATIENT COMFORT
- ▶ EASY TO USE ABS MOULDED ARM REST
- GLASS SPITTOON WITH AUTO CLOSE WATER CUP MER AND SPITTOON FLUSH
- ▶ 4 LED LAMP WITH SENSOR HAVING ILLUMINATION BETWEEN 7000-32000
- ▶ MULTIPURPOSE UTILITY TRAY FOR SCALER, GLASS & NAPKIN HOLDER
- LARGE INSTRUMENT TRAY WITH STAINLESS STEEL TOP & PREUMATOC LOCKING SYSTEM
- ▶ IN DELIVERY UNIT: PROVISION FOR 2 ANSOTOX, 1 MICROMOTER &1 THREEWAY SYARIGE
- ▶ IN ASSISTANT UNIT: 1 THREEWAY SYRINGE & A LOW SUCTION
- X-RAY VIEWER WITH WHITE LIGHT
- ▶ MULTIFUNTIONAL FOOT CONTROL AND DOUBLE ARTICULATOR HEADREST FOR HANDICAPPED & PAEDIATRIC PATIENTS

Accompanying the Recovery Journey

KEY FEATURE

THE DENTAL CHAIR AQUA IS AN ADVANCED AND ERGONOMICALLY DESIGNED TREATMENT UNIT TAILORED FOR MODERN DENTAL CLINICS. FEATURING ELECTRICAL OPERATION WITH TWO PRE-SET POSITIONS—ZERO AND SUPINE, IT ENSURES EFFORTLESS PATIENT POSITIONING FOR OPTIMAL TREATMENT EFFICIENCY. DESIGNED WITH SEAMLESS AND SOFT UPHOLSTERY, IT ENHANCES PATIENT COMFORT DURING PROLONGED PROCEDURES, WHILE THE ABS-MOULDED ARMREST ADDS CONVENIENCE AND SUPPORT. THE CHAIR IS EQUIPPED WITH A GLASS SPITTOON FEATURING AN AUTO-CLOSE WATER CUP FILLER AND SPITTOON FLUSH, ENSURING SUPERIOR HYGIENE STANDARDS. A 4-LED SENSOR LAMP WITH ADJUSTABLE ILLUMINATION BETWEEN 7,000-32,000 LUX PROVIDES CLEAR VISIBILITY, AIDING IN PRECISION TREATMENTS. THE MULTIPURPOSE UTILITY TRAY IS THOUGHTFULLY DESIGNED TO HOLD ESSENTIAL TOOLS LIKE A SCALER, GLASS, AND NAPKIN, ENSURING EASY ACCESSIBILITY. A LARGE STAINLESS STEEL INSTRUMENT TRAY WITH A PNEUMATIC LOCKING SYSTEM FURTHER ENHANCES FUNCTIONALITY. THE DELIVERY UNIT IS EQUIPPED WITH PROVISIONS FOR TWO AIROTORS, ONE MICROMOTOR, AND ONE THREE-WAY SYRINGE, WHILE THE ASSISTANT UNIT INCLUDES A THREE-WAY SYRINGE AND LOW SUCTION FOR SEAMLESS DENTAL PROCEDURES. ADDITIONALLY, THE BUILT-IN X-RAY VIEWER WITH WHITE LIGHT ASSISTS IN PRECISE DIAGNOSIS. A MULTIFUNCTIONAL FOOT CONTROL AND DOUBLE-ARTICULATED HEADREST MAKE IT SUITABLE FOR PEDIATRIC AND HANDICAPPED PATIENTS, ENSURING MAXIMUM ADAPTABILITY. WITH ITS INNOVATIVE FEATURES, USER-FRIENDLY CONTROLS, AND ERGONOMIC DESIGN, THE DENTAL CHAIR AQUA IS A PERFECT ADDITION TO ANY MODERN DENTAL PRACTICE, ENHANCING BOTH EFFICIENCY AND PATIENT CARE.

Create a Safer Environment



Procare

SEC NO. 7 , PLOT NO. 146A, PCNTDA MIDC
BHOSARI, PUNE - 411039, MAHARASHTRA
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare