



DENTAL CHAIR PLATINUM

www.procaremedfurn.com

ELECTRICALLY OPERABLE DENTAL TREATMENT UNIT WITH 5 PROGRAMS: 3 PROGRAMMABLE & 2 FIXED: ZERO & GARGLE POSITION


- ▶ SEAMLESS & SOFT UPHOLSTERY FOR PATIENT COMFORT
- ▶ EASY TO USE ABS MOULDED ARM REST
- ▶ MOVABLE GLASS SPITTOON WITH AUTO CLOSE WATER CUP FILLER
- ▶ 6 LED LAMP WITH SENSOR HAVING ILLUMINATION BETWEEN 7000-32000 LUX
- ▶ MULTIPURPOSE UTILITY TRAY FOR SCALER, GLASS & NAPKIN HOLDER
- ▶ FEATHER TOUCH CHAIR CONTROL PRESENT AT BOTH SIDE: DELIVERY UNIT & IN ASSISTANT SIDE
- ▶ LARGE INSTRUMENT TRAY WITH STAINLESS STEEL TOP & PNEUMATIC LOCKING SYSTEM
- ▶ IN DELIVERY UNIT: PROVISION FOR 3 AIROTOR, 1 MICROMOTOR & 1 THREWAY SYNINGE
- ▶ IN ASSISTANT UNIT: 1 THREWAY SYRINGE & A LOW SUCTION
- ▶ X-RAY VIEWER WITH WHITE LIGHT
- ▶ MULTIFUNTIONAL FOOT CONTROL AND DOUBLE ARTICULATED HEADREST FOR HANDICAPPED & PAEDIATRICPATIENTS

Accompanying the Recovery Journey

KEY FEATURE

THE DENTAL CHAIR PLATINUM IS A STATE-OF-THE-ART SOLUTION DESIGNED TO ENHANCE PATIENT COMFORT AND STREAMLINE DENTAL PROCEDURES. FEATURING ELECTRICAL OPERATION WITH FIVE PRE-SET PROGRAMS—THREE PROGRAMMABLE AND TWO FIXED (ZERO & GARGLE POSITIONS)—THIS CHAIR ENSURES SEAMLESS ADJUSTMENTS FOR OPTIMAL PATIENT POSITIONING. ITS SOFT AND SEAMLESS UPHOLSTERY PROVIDES SUPERIOR COMFORT, WHILE THE ABS-MOULDED ARMREST ADDS EASE OF USE. THE CHAIR IS EQUIPPED WITH A MOVABLE GLASS SPITTOON WITH AN AUTO-CLOSE WATER CUP FILLER, PROMOTING HYGIENE AND CONVENIENCE. A 6-LED SENSOR LAMP WITH ADJUSTABLE ILLUMINATION BETWEEN 7,000-32,000 LUX ENSURES PRECISION DURING TREATMENTS. THE MULTIPURPOSE UTILITY TRAY EFFICIENTLY HOLDS ESSENTIAL TOOLS LIKE A SCALER, GLASS, AND NAPKIN FOR EASY ACCESS. FOR ENHANCED USABILITY, THE CHAIR FEATURES FEATHER-TOUCH CONTROLS ON BOTH THE DELIVERY UNIT AND ASSISTANT SIDE, ALONG WITH A LARGE STAINLESS STEEL INSTRUMENT TRAY WITH A PNEUMATIC LOCKING SYSTEM FOR ORGANIZED TOOL MANAGEMENT. THE DELIVERY UNIT INCLUDES PROVISIONS FOR THREE AIROTORS, ONE MICROMOTOR, AND ONE THREE-WAY SYRINGE, WHILE THE ASSISTANT UNIT IS EQUIPPED WITH A THREE-WAY SYRINGE AND LOW SUCTION. ADDITIONALLY, THE CHAIR BOASTS A BUILT-IN X-RAY VIEWER WITH WHITE LIGHT, AIDING IN CLEAR DIAGNOSTICS. A MULTIFUNCTIONAL FOOT CONTROL AND A DOUBLE-ARTICULATED HEADREST ENSURE COMFORT AND ACCESSIBILITY FOR PEDIATRIC AND HANDICAPPED PATIENTS. COMBINING CUTTING-EDGE FEATURES, ERGONOMIC DESIGN, AND SEAMLESS OPERATION, THE DENTAL CHAIR PLATINUM IS THE PERFECT CHOICE FOR MODERN DENTAL CLINICS SEEKING EFFICIENCY, PRECISION, AND ENHANCED PATIENT CARE.

Create a Safer Environment



Procare

SEC NO. 7 , PLOT NO. 146A, PCNTDA MIDC
BHOSARI, PUNE - 411039, MAHARASHTRA
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.