

DERMA CHAIR ELECTRIC

www.procaremedfurn.com

GENERAL USE :-MEDICAL BED

MATERIAL :- SYNTHETIC LEATHER PVC/PU LEATHER, SPONGE

SIZE :-70LX22WX26-36H INCH

LOADING :- 200KGS

REMOTE :- PROVIDED

MOTORS :- 3 MOTORS

PILLOW :- SEPARATE

COLOR :-VARIOUS COLOR AVAILABLE

USAGE :- HAIR TRANSPLANT, SKIN COSMETOLOGY DERMA

CLINIC & HOSPITAL

OPTIONAL :- BATTERY BACKUP

Accompanying the Recovery Journey

KEY FEATURE

ENHANCE PATIENT EXPERIENCE AND PROCEDURAL EFFICIENCY WITH OUR ELECTRIC DERMA CHAIR, DESIGNED FOR HAIR TRANSPLANT, SKIN COSMETOLOGY, DERMA CLINICS, AND HOSPITALS. ENGINEERED FOR MAXIMUM COMFORT, DURABILITY, AND ADJUSTABILITY, THIS CHAIR IS IDEAL FOR PRECISION TREATMENTS.

KEY FEATURES:

- 1.ERGONOMIC & SPACIOUS DESIGN: WITH DIMENSIONS 70L X 22W X 26-36H INCHES, IT PROVIDES AMPLE SPACE FOR PATIENT COMFORT.
 2. HIGH-QUALITY MATERIAL: CRAFTED FROM SYNTHETIC LEATHER (PVC/PU) AND SPONGE PADDING, ENSURING A SOFT, EASY-TO-CLEAN SURFACE.
 3. ROBUST & DURABLE: SUPPORTS UP TO 200 KG, OFFERING SUPERIOR STABILITY FOR VARIOUS TREATMENTS.
 4. ADVANCED MOTORIZED ADJUSTMENTS: 3 POWERFUL MOTORS FOR SEAMLESS HEIGHT, RECLINE, AND TILT ADJUSTMENTS.
 5. USER-FRIENDLY REMOTE CONTROL: ENABLES EFFORTLESS POSITIONING ADJUSTMENTS FOR BOTH PATIENTS AND PRACTITIONERS.
 6. SEPARATE PILLOW: ENHANCES HEAD AND NECK SUPPORT FOR EXTENDED PROCEDURES.
 7. CUSTOMIZABLE AESTHETICS: AVAILABLE IN VARIOUS COLORS TO SUIT CLINIC OR HOSPITAL INTERIORS.
- OPTIONAL BATTERY BACKUP: ENSURES UNINTERRUPTED OPERATION IN CASE OF POWER OUTAGES.
- IDEAL FOR COSMETIC PROCEDURES, DERMATOLOGY TREATMENTS, AND HAIR TRANSPLANTS, THIS ELECTRIC DERMA CHAIR DELIVERS EFFICIENCY, COMFORT, AND MODERN TECHNOLOGY, MAKING IT A MUST-HAVE FOR PROFESSIONAL CLINICS.
- UPGRADE YOUR DERMA CLINIC OR HOSPITAL WITH THE BEST ELECTRIC DERMA CHAIR TODAY!


Create a Safer Environment



CE

Procare

SEC NO. 7 , PLOT NO. 146A, PCNTDA MIDC
BHOSARI, PUNE - 411039, MAHARASHTRA
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare